

JETT PLASMA For Her II

FOR GYNAECOLOGY



REJUVENATES
THE WORLD!

Mechanism

Possitive effect of JETT PLASMA For Her II is caused by combination of heat, membrane depolarization and low threshold reversible electroporation. The passage of direct current generates heat, which activates the formation of new collagen connective tissue, leading to a rejuvenation and stronger and more flexible tissue structure. In addition to the heat effect, unique direct current technology causes membrane depolarization and reversible electroporation which reduces the barrier and enable easier passage of ions and molecules through cell membranes.

Advantages

- Safe, comfortable, painless
- Short treatment time
- Non-surgical
- No need for anesthesia
- No hormones
- No scars
- 3 treatments with an interval of 10-14 days

JETT PLASMA For Her II
with Plasma Pen III/G



UNIQUE APPLICATORS DEVELOPED FOR PATIENT COMFORT DELIVERED IN A STERILE STATE

Plasma Pen III/G
with **Flat Applicator 10 mm**
- for treatment of vulvar laxity



G-Applicator Wave
- for treatment of vulvovaginal laxity



G-Applicator 30 mm
- for treatment of vaginal laxity



JettPlasmaDevices

COMPEX, spol. s r.o., Palackého tř. 924/105, 612 00 Brno, CZ | info@jett.eu | www.jett.eu